

Analysis Date:

Re-test Date:

ABEMACICLIB D10

Identification


Chemical Name : N-(5-((4-Ethylpiperazin-1-yl)methyl)pyridin-2-yl)-5-fluoro-4-(4-fluoro-2-(methyl-d3)-1-(propan-2-yl-d7)-1H-benzo[d]imidazol-6-yl)pyrimidin-2-amine

CAT No : ALL-A07526

CAS No : 2088650-71-7

Molecular Formula : C₂₇H₂₂D₁₀F₂N₈

Molecular Weight : 516.7



Analytical Information

Batch Code : ALL-A07526

HPLC Purity : 98.00 %

Solubility : MeOH: ACN

Potency : N.A

Appearance of Product : Off White Solid

Mass : Confirm

Long Term Storage : 2-8°C

IR Analysis : Confirm

Weight Loss By TGA : N.A

¹H NMR : Confirm

Residue Of Ignition : N.A

Additional Information

%Potency = [100 -(Inorganic Impurities % + Water %) × Chromatographic Purity%]/100 = N.A

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Swapnil Kausatkar Jr. Executive	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, ¹H NMR, IR

Shipping Condition : All Product are stable to be shipped at room temperature, unless otherwise specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E), Mumbai, Maharashtra 421203