

Analysis Date:

Re-test Date:

N-ACETYL-D-CYSTEINE

Identification


Chemical Name : Acetyl-D-cysteine

CAT No : ALL-A07743

CAS No : 26117-28-2

Molecular Formula : C₅H₉NO₃S

Molecular Weight : 163.2



Analytical Information

Batch Code : ALL-A07743

HPLC Purity : 98.00 %

Solubility : MeOH: ACN

Potency : N.A

Appearance of Product : Off White Solid

Mass : Confirm

Long Term Storage : 2-8⁰C

IR Analysis : Confirm

Weight Loss By TGA : N.A

¹HNMR : Confirm

Residue Of Ignition : N.A

Additional Information

%Potency = [100 -(Inorganic Impurities % + Water %) × Chromatographic Purity%]/100 = N.A

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Swapnil Kausatkar Jr. Executive	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, 1H NMR, IR

Shipping Condition : All Product are stable to be shipped at room temperature, unless otherwise specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E), Mumbai, Maharashtra 421203