
Certificate Of Analysis

Analysis Date: 10-05-2025 **Re-test Date:** 10-05-2028

AFATINIB IMPURITY 19

Identification

Chemical Name : (E)-N-(4-((3-Chloro-4-fluorophenyl)amino)-7-methoxyquinazolin-6-yl)-4-(dimethylamino)but-2-enamide



CAT No : ALL-A07975

CAS No : 1245555-42-3

Molecular Formula : C21H21ClFN5O2

Molecular Weight : 429.9

Analytical Information

Batch Code	: ALL-A07975	**HPLC Purity**	: 98.06%
Solubility	: MeOH: ACN	**Potency**	: 96.82%
Appearance of Product	: Off White Solid	**Mass**	: Confirm
Long Term Storage	: 2-8 ^{0}C	**IR Analysis**	: Confirm
Weight Content By TGA	: 0.583%	**1HNMR**	: Confirm
Residue Of Ignition	: 0.673%		

Additional Information

%Potency = [100 - (Inorganic Impurities % + Water %) × Chromatographic Purity%]/100 =

[100 – (0.583+0.673) × 98.06]/100 = 96.82%

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Vipul Khadse Jr. Executive	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, 1H NMR, IR, TGA

Shipping Condition : All Product are stable to be shipped at room temperature, unless otherwise specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E),Mumbai, Maharashtra 421203