Certificate Of Analysis

Analysis Date: 13-05-2025 **Re-test Date:** 13-05-2028

3-(4-METHOXYBENZOYL) PROPIONIC ACID

Identification


Chemical Name : 4-(4-methoxyphenyl)-4-oxobutanoic acid

CAT No : ALL-A08099

CAS No : 3153-44-4

Molecular Formula : C11H12O4

Molecular Weight : 208.2



Analytical Information

Batch Code	: ALL-A08099	**HPLC Purity**	: 98.02%
Solubility	: MeOH: ACN	**Potency**	: 96.83%
Appearance of Product	: Off White Solid	**Mass**	: Confirm
Long Term Storage	: 2-8 ^{0}C	**IR Analysis**	: Confirm
Weight Content By TGA	: 0.523%	**1HNMR**	: Confirm
Residue Of Ignition	: 0.681%		

Additional Information

%Potency = [100 - (Inorganic Impurities % + Water %) × Chromatographic Purity%]/100 =

[100 – (0.523+0.681× 98.02]/100 = 96.83%

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Vipul Khadse Jr. Executive	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, 1H NMR, IR, TGA

Shipping Condition : All Product are stable to be shipped at room temperature, unless otherwise specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E),Mumbai, Maharashtra 421203