Certificate Of Analysis

Analysis Date: 13-05-2025 **Re-test Date:** 13-05-2028

AGOMELATINE 7-DES METHYL IMPURITY

Identification


Chemical Name : (2-(7-Hydroxynaphthalen-1-yl)ethyl)acetamide

CAT No : ALL-A04785

CAS No : 52302-45-9

Molecular Formula : C14H15NO2

Molecular Weight : 229.27



Analytical Information

Batch Code	: ALL-A04785	**HPLC Purity**	: 98.04%
Solubility	: MeOH: ACN	**Potency**	: 96.65%
Appearance of Product	: Off White Solid	**Mass**	: Confirm
Long Term Storage	: 2-8 ^{0}C	**IR Analysis**	: Confirm
Weight Content By TGA	: 0.581%	**1HNMR**	: Confirm
Residue Of Ignition	: 0.834%		

Additional Information

%Potency = [100 - (Inorganic Impurities % + Water %) × Chromatographic Purity%]/100 =

[100 – (0.581+0.834× 98.04]/100 = 96.65%

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Vipul Khadse Jr. Executive	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, 1H NMR, IR, TGA

Shipping Condition : All Product are stable to be shipped at room temperature, unless otherwise specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E),Mumbai, Maharashtra 421203