Certificate Of Analysis

Analysis Date: 13-05-2025 **Re-test Date:** 13-05-2028

AGOMELATINE IMPURITY 16

Identification


Chemical Name : 2-(7-Methoxynaphthalen-1-yl)acetic acid

CAT No : ALL-A08106

CAS No : 6836-22-2

Molecular Formula : C13H12O3

Molecular Weight : 216.2



Analytical Information

Batch Code	: ALL-A08106	**HPLC Purity**	: 98.02%
Solubility	: MeOH: ACN	**Potency**	: 96.94%
Appearance of Product	: Off White Solid	**Mass**	: Confirm
Long Term Storage	: 2-8 ^{0}C	**IR Analysis**	: Confirm
Weight Content By TGA	: 0.689%	**1HNMR**	: Confirm
Residue Of Ignition	: 0.410%		

Additional Information

%Potency = [100 - (Inorganic Impurities % + Water %) × Chromatographic Purity%]/100 =

[100 – (0.689+0.410× 98.02]/100 = 96.94%

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Vipul Khadse Jr. Executive	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, 1H NMR, IR, TGA

Shipping Condition : All Product are stable to be shipped at room temperature, unless otherwise specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E),Mumbai, Maharashtra 421203