Certificate Of Analysis

Analysis Date: 15-05-2025 **Re-test Date:** 15-05-2028

AGOMELATINE UREA

Identification


Chemical Name : N-(2-(7-methoxynaphthalen-1- yl)ethyl)acetamide urea

CAT No : ALL-A04788

CAS No : 1366181-82-9

Molecular Formula : C16H21N3O3 (Free Base)

Molecular Weight : 303.36 (Free Base)



Analytical Information

Batch Code	: ALL-A04788	**HPLC Purity**	: 98.01%
Solubility	: MeOH: ACN	**Potency**	: 96.66%
Appearance of Product	: Off White Solid	**Mass**	: Confirm
Long Term Storage	: 2-8 ^{0}C	**IR Analysis**	: Confirm
Weight Content By TGA	: 0.594%	**1HNMR**	: Confirm
Residue Of Ignition	: 0.774%		

Additional Information

%Potency = [100 - (Inorganic Impurities % + Water %) × Chromatographic Purity%]/100 =

[100 – (0.594+0.774× 98.01]/100 = 96.66%

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Vipul Khadse Jr. Executive	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, 1H NMR, IR, TGA

Shipping Condition : All Product are stable to be shipped at room temperature, unless otherwise specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E),Mumbai, Maharashtra 421203