Certificate Of Analysis

Analysis Date: 20-05-2025 **Re-test Date:** 20-05-2028

ALBENDAZOLE EP IMPURITY I

Identification


Chemical Name : Methyl N-(5-Propoxy-1H-benzimidazol-2-yl)carbamate

CAT No : ALL-ABZ-73

CAS No : 20559-55-1

Molecular Formula : C12H15N3O3

Molecular Weight : 249.27



Analytical Information

Batch Code	: ALL-ABZ-73	**HPLC Purity**	: 98.09%
Solubility	: MeOH: ACN	**Potency**	: 97.03%
Appearance of Product	: Off White Solid	**Mass**	: Confirm
Long Term Storage	: 2-8 ^{0}C	**IR Analysis**	: Confirm
Weight Content By TGA	: 0.478%	**1HNMR**	: Confirm
Residue Of Ignition	: 0.593%		

Additional Information

%Potency = [100 - (Inorganic Impurities % + Water %) × Chromatographic Purity%]/100 =

[100 – (0.478+0.593× 98.09]/100 = 97.03%

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Vipul Khadse Jr. Executive	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, 1H NMR, IR, TGA

Shipping Condition : All Product are stable to be shipped at room temperature, unless otherwise specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E),Mumbai, Maharashtra 421203