Certificate Of Analysis

Analysis Date: 24-12-2025 **Re-test Date:** 24-12-2028

SERTRALINE-D6 HYDROCHLORIDE

Identification


Chemical Name : SERTRALINE-D6

CAT No : ALL-D01155

CAS No : **N. A.**

Molecular Formula : C17H11D6Cl2N • HCl

Molecular Weight : 348.73



Analytical Information

Batch No.	: ALL-D01155	**HPLC Purity**	: 98.00%
Solubility	: USP Diluent / EP Diluent (MEOH)	**% D Atom**	: 99.00 %
Appearance of Product	: Off White Solid	**Potency**	: 98.00%
Weight Loss By TGA	: 0.115%	**Mass**	: Confirm
Residue Of Ignition	: 0.148%	**IR Analysis**	: Confirm
Long Term Storage	: -20°C	**1HNMR**	: Confirm

Additional Information

%Potency = [100 - (Weight Loss By TGA % + Residue Of Ignition %) × Chromatographic Purity%]/100 =

[100 – (:0.115+0.148) × 98.00]/100 = 97.74%

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr.Vipul khadase Jr. Executive	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, 1H NMR, IR, TGA

Shipping Condition : All Product are stable to be shipped at room temperature, unless specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E),Mumbai, Maharashtra 421203